AT04511

MEDYCZNY PAS ROZGRZEWAJĄCY INSTRUKCJA UŻYTKOWANIA
MEDICAL WARMING BELT INSTRUCTIONS FOR USE
GEBRAUCHSANWEISUNG FÜR MEDIZINISCHEN WÄRMEGURT
NÁVOD K POUŽITÍ ZDRAVOTNÍHO ZAHŘÍVACÍHO PÁSU
NÁVOD NA POUŽITIE ZDRAVOTNÉHO ZAHREJOVACIEHO PÁSU
NOTICE D'UTILISATION DE LA CEINTURE CHAUFFANTE MEDICALE
MEDISCHE VERWARMINGSRIEM INSTRUCTIES VOOR GEBRUIK
INSTRUCCIONES DE USO DEL CINTURÓN DE CALENTAMIENTO MÉDICO
FASCIA RISCALDANTE MEDICA ISTRUZIONI PER L'USO
ANVÄNDNINGSSINSTRUKTIONER FÖR MEDICINSK VÄRMEBÄL

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INTENDED USE AND PERFORMANCE CHARACTERISTICS OF THE PRODUCT

The medical warming belt filled with cherry seeds is intended for people complaining of pain and tension in the lumbar region of the spine, it helps to relieve pain and muscle tension. The product can be used hot or cold.

Waist circumference: 80-125 cm

CONTRAINDICATIONS FOR USE OF THE PRODUCT

Do not use if you are allergic to any of the materials included in the product (cotton/cherry seeds).

CAUTION

In the event of a product-related „serious incident“ that directly or indirectly led to, may have led to, or may lead to, any of the following:

- death of a patient, user or other person or
- temporary or permanent deterioration of the patient's, user's or other person's health or
- a serious threat to public health

this „serious incident“ should be reported to the manufacturer and to the competent authority of the Member State where the user or patient resides. In the case of Poland, the competent authority is the Office for Registration of Medicinal Products, Medical Devices and Biocidal Products.

CAUTION

In the event of pain, allergic reactions or other disturbing, unclear for the user symptoms related to the use of the medical device, consult a healthcare professional.

**Warnings and Precautions:**

Contains small parts which present a choking hazard to children.
Prohibited use by children under 3 years of age.

ATTENTION: it is forbidden to use the product in a way other than for its intended purpose.

NOTE: Before use, check the product for damage, such as holes in the material.

PRODUCT PREPARATION

For heat: Put the product in the microwave oven, oven, put on the radiator or any source of heat.

Cold: place the product in the freezer for about 30 minutes

CLEANING

In case of dirt, the cover should be washed by hand.

STORAGE CONDITIONS: Store in a dry place, away from direct sunlight and moisture.

PACKAGE INCLUDES: 1 cover, 4 inserts.

TARGET GROUP OF PATIENTS

The health care professional indicates the use of the device to adults and children on his own responsibility, taking into account the available variants/sizes/necessary functions/size and indications, taking into account the information provided by the manufacturer.

METHOD OF DISPOSAL OF THE PRODUCT AFTER WITHDRAWAL FROM USE

After the device is decommissioned, the medical device can be disposed of as ordinary municipal waste.

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WARRANTY CARD

Model:
Serial number:

Point of sale stamp
and a legible signature of the Seller:

1. ANTAR Sp. J., 03-068 Warsaw, ul. Zawisłarska 43, hereby grants a 12-month warranty on the rollator, starting from the date of purchase by the Buyer. Products financed by the National Health Fund are covered with the extended warranty period, which covers half of the useful life specified in the Regulation on Guaranteed Services, provided that the useful life is specified in the Regulation in months or years.
2. During the warranty period ANTAR undertakes to carry out the repairs necessary to restore proper functioning of the product, within 14 days from the date of receipt of the product by the service centre, together with the proof of purchase and the warranty card.
3. ANTAR reserves the right to replace the product if the repair is deemed uneconomic. Repair or replacement of faulty products does not extend the warranty period.
4. If the product is found to be malfunctioning, send it immediately to ANTAR (as specified in point 1) or contact the point of sale.
5. Before use, read the instructions carefully and follow them. If the product is used improperly and contrary to recommendations, the warranty will be irrevocably lost.
6. The warranty does not cover damage caused by external mechanical force, load exceeding the maximum permissible weight of the user as described in the manual, as well as contact of the product with alcohol, greasy substances or petrol. The warranty does not cover parts subject to normal wear and tear during use, such as the bearing, rubber parts and others.
7. The product delivered to the service in connection with the complaint cannot cause an epidemic.
8. THIS WARRANTY IS ONLY VALID TOGETHER WITH THE PROOF OF PURCHASE (BILL, RECEIPT, INVOICE). IN CASE OF A COMPLAINT, A PROOF OF PURCHASE MUST BE ATTACHED. IN ORDER TO EXTEND WARRANTY PERIOD FOR PRODUCTS FINANCED BY THE NATIONAL HEALTH FUND, ATTACH A COPY OF THE ORDER FOR THE SUPPLY OF MEDICAL DEVICES BEING ORTHOPAEDIC ITEMS AND AIDS.
9. Complaints submitted without the proof of purchase and a warranty card with the product serial number will not be considered.

 **ANTAR Sp. J.**
03-068 Warsaw, ul. Zawisłarska 43
Tel. 22 518 36 00, Fax 22 518 36 30
www.antar.net

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