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ENG

GENERAL WARNING

Thank you for purchasing our product, please read the instructions for use carefully. Do not use the product before reading and understanding this manual. If you do not understand the warnings, cautions and recommendations, contact your healthcare professional or retailer to avoid damaging the product.

APPLICATION/INDICATIONS


Prevention of swelling and stabilization in the case of: soft tissue injuries in the area of the ankle joint, foot and distal part of the shin, sprains of the ankle joint, metatarsal and tarsus of the second and third degree, treatment of stable fractures in the area of the ankle joint and metatarsal, after surgical procedures in the above-mentioned area, changes inflammation of the Achilles tendon and foot extensor and flexor tendons.

CONTRAINDICATIONS

physical or mental limitations (e.g. visual impairment) that prevent safe handling of the product.


USE

1. Place your foot and leg into the foam insole and tighten so that the insole is snug against your foot.
2. Use cushions to remove empty space in the insert or between the insert and the housing.
3. Place the foot wrapped in the insole in the shoe and align it to the side axis on both sides. Remove the plastic caps from the supports to slide them into the foam insert.
4. Tighten all straps one by one and fasten them.


CLOCK

1. Turn the plastic cap to adjust the clock.
2. Move the clock to the desired setting.
3. Lock in the right position with the metal lock.
4. Turn the plastic cap back to secure the clock

ADJUSTING THE PRESSURE IN THE AIR CHAMBER

1. Turn the button to the "inflate" position, then press the pump to inflate the chamber.
2. Turn the button to the "deflate" position to deflate the chamber.

SIZES AT53003, AT53004

shoe size:
<35.5 - XS
36-39 - p
40-44 - M
44-46.5 - L
>46.5 - XL

SCOPE OF DELIVERY

Orthosis, screw to adjust the clock, spare locks to adjust the flexion angle.

CLEANING AND MAINTENANCE

It is recommended to wash the product by hand with a mild detergent in cool water. Do not machine wash. Do not use strong detergents, bleach or abrasive objects for cleaning. Dry flat in the shade, do not expose to direct sunlight. Do not tumble dry. Do not dry clean. Do not iron. Remove all plastic and steel parts before washing.

STORAGE

The product should be stored in a dry, cool room, away from moisture and direct sunlight.

ATTENTION:

In the event of a product-related "serious incident" that directly or indirectly led to, may have led to, or may lead to, any of the following:

- a) death of a patient, user or other person or
- b) temporary or permanent deterioration of the patient's, user's or other person's health or
- (c) a serious threat to public health

this "serious incident" should be reported to the manufacturer and to the competent authority of the Member State where the user or patient resides. In the case of Poland, the competent authority is the Office for Registration of Medicinal Products, Medical Devices and Biocidal Products.

ATTENTION:

In the event of pain, allergic reactions or other disturbing, unclear for the user symptoms related to the use of the medical device, consult a healthcare professional.

ATTENTION: it is forbidden to use the product in a way other than for its intended purpose

NOTE: The manufacturer is not responsible for damage caused by neglecting the rules of cleaning and maintenance, or resulting from failure to comply with the recommendations contained in this manual.

TARGET PATIENT GROUP

Persons suffering from diseases, dysfunctions or injuries for the treatment/rehabilitation or compensation of which this device is intended (see section on intended use of the device in this instruction manual). The device can be purchased by the user on their own or on the recommendation of a doctor, therapist or other specialist. Whether you purchase the device yourself or on the recommendation of a doctor/therapist/other specialist, you must take into account the available sizes/necessary functions and variants of the device, the indications and contraindications for use and the information provided by the manufacturer.

HOW TO DISPOSAL OF THE PRODUCT

After the device is decommissioned, the medical device can be disposed of as ordinary municipal waste

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Model:
Serial number:

Point of sale stamp

and a legible signature of the Seller:

1. ANтар Sp. J., 03-068 Warsaw, ul. Zawiślаńska 43, hereby grants a 12-month warranty on the product, starting from the date of purchase by the Buyer.
2. During the warranty period ANтар undertakes to carry out the repairs necessary to restore proper functioning of the product, within 14 days from the date of receipt of the product by the service centre, together with the proof of purchase and the warranty card.
3. ANтар reserves the right to replace the product if the repair is deemed uneconomic. Repair or replacement of faulty products does not extend the warranty period.
4. If the product is found to be malfunctioning, please contact your retailer/distributor or (if it is impossible) send it to ANтар (as specified in point 1).
5. Before use, read the instructions carefully and follow them. If the product is used improperly and contrary to recommendations, the warranty will be irrevocably lost.
6. The warranty does not cover damage caused by external mechanical force, load exceeding the maximum allowed weight of the user, as described in the manual, as well as contact of the product with alcohol, greasy substances or petrol. The warranty does not cover parts which are subject to normal wear or tear during use, such as the bearing, rubber parts and others.
7. The product delivered to the service in connection with the complaint cannot cause an epidemic hazard.
8. THIS WARRANTY IS ONLY VALID TOGETHER WITH THE PROOF OF PURCHASE (BILL, RECEIPT, INVOICE). IN CASE OF A COMPLAINT, A PROOF OF PURCHASE MUST BE ATTACHED.
9. Complaints submitted without the proof of purchase and a warranty card with the product serial number will not be considered.

ANTAR Sp. J.
03-068 Warsaw, ul. Zawiślаńska 43
Tel. 22 518 36 00, Fax 22 518 36 30
www.antar.net

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