AT53005, AT53006, AT53008, AT53030
ORTEZA SZTYWNA Z TWORZYWA SZTUCZNEGO NA GOLEŃ I STOPE
INSTRUKCJA UŻYTKOWANIA
PLASTIC RIGID SHIN AND FOOT ORTHOSIS USER MANUAL
GEBRAUCHSANWEISUNG FÜR STARRE SCHIENBEIN- UND FUSSORTHESEN
AUS KUNSTSTOFF
PLASTOVÁ PEVNÁ ORTÉZA NA HOLEŇ A CHODIDLO UŽIVATELSKÁ PŘÍRUČKA
PLASTOVÁ PEVNÁ ORTÉZA NA HOLEŇ A CHODIDLO POUŽÍVATEĽSKÁ
PRÍRUČKA
ORTHÈSE RIGIDE EN PLASTIQUE POUR LA JAMBE ET LE PIED
MODE D'EMPLOI
KUNSTSTOF RIGIDE ORTHESE VOOR ONDERBEEN EN VOET
GEBRUIKSAANWIJZING
ORTESIS RÍGIDA DE PLÁSTICO PARA LA PIERNA Y EL PIE
MANUAL DE INSTRUCCIONES
ORTESI RIGIDA IN PLASTICA PER GAMBA E PIEDE
MANUALE D'USO
STYV PLASTORTOS FÖR UNDERBEN OCH FOT
BRUKSANVISNING

MD

PL

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ENG

Thank you for purchasing our product. Please carefully read the user manual. Do not use the product before reviewing and understanding this manual. If you do not understand the warnings, notes, and recommendations, please consult a healthcare professional or the seller to avoid product damage or personal injury.

NOTE: Check all parts for damage that may have occurred during transport. If such damage is noticed, DO NOT use the product. More information is available from the seller.

PURPOSE The AT53005/AT53006 orthosis is a high-quality single-chamber rigid brace made of plastic, designed as an essential tool for treating various injuries and conditions related to the ankle joint, foot, and shin. Carefully crafted, this orthosis provides excellent stabilization and support for injured tissues, thereby effectively supporting the recovery process. It is lined with soft-touch fabric, lightweight, and ergonomic, guaranteeing comfortable wear over long periods. Additionally, pneumatic cushions allow an ideal fit to the shape of the foot and shin, ensuring comfort and excellent immobilization during therapy. Applications:

Treatment of bone fractures near the ankle and foot

Soft tissue injuries, including the Achilles tendon

Ankle sprains

The AT53008/AT53030 is a single-chamber rigid orthosis made of plastic, dedicated to individuals dealing with various injuries and conditions of the ankle joint and foot. The short version of the orthosis provides excellent stabilization and support for the foot, ankle, and lower leg. Its features make it an ideal solution for dislocations, sprains, fractures, and post-operative treatment. Comfortable fastening, soft-touch fabric, and lightweight design make it an effective yet comfortable therapeutic tool. Additionally, a solid, shock-absorbing sole and pneumatic cushions guarantee an excellent fit and user comfort, serving as an alternative to traditional plaster casts. Applications:

Ankle dislocations and sprains

Post-operative treatment

Fractures near the ankle joint

Foot fractures

Metatarsal bone fractures

Achilles tendon rupture and inflammation

Alternative to standard plaster casts

CONTRAINDICATIONS Physical or mental limitations (e.g., visual impairment) that prevent safe handling of the product.

Skin infections and open wounds – Orthoses should not be used on areas with wounds, abrasions, skin infections, ulcers, or other dermatological problems. Lack of air access and pressure may aggravate the infection.

Swelling and venous insufficiency – In cases of swelling or circulation issues, using a rigid orthosis may cause discomfort and worsen inflammation.

Limb ischemia – People with ischemic problems in the lower limbs, such as atherosclerosis, may experience increased discomfort with rigid orthoses as they may affect circulation.

Osteoporosis and bone fragility – In individuals with advanced osteoporosis, using a rigid orthosis may lead to compression fractures or other injuries due to bone pressure.

Acute injury phase with severe pain and inflammation – The orthosis may not be recommended if the injury is in an acute phase with significant swelling and pain, as it may cause additional pressure and exacerbate symptoms.

Contractures or joint deformities – Individuals with advanced deformities or contractures may not be able to fit into a rigid orthosis, which limits its functionality and effectiveness.

Neurological problems and sensory disturbances – Patients with sensory disturbances (e.g., neuropathy) may not feel if the orthosis is too tight or causes other discomfort, increasing the risk of complications.

TARGET PATIENT GROUP People suffering from diseases, dysfunctions, or injuries for the treatment/rehabilitation or compensation of which this device is intended (see the product's intended use in this manual). The product can be purchased by the user independently or on the recommendation of a doctor, therapist, or another specialist. Whether the product is purchased independently or recommended by a doctor/therapist/other specialist, consider available sizes/necessary functions and variants, indications and contraindications for use, as well as information provided by the manufacturer.

USE





1 Place the foot and leg in the foam insole and tighten so that the insole is snug against the foot.

2. for a better fit, place cushions inside the insole or between the insole and the casing.

Place the foot wrapped in the insole in the shoe and pull the vertical supports on both sides slightly outwards.


outwards. The lateral axis of the shin should be aligned on both sides to the supports. Then release the supports so that they stick to the foam insole. Slide the first strap from the top into the holes on both sides of the housing.

4. tighten all hooks and straps in order from the front of the foot to the top of the leg.

ADJUSTING AIR CHAMBER PRESSURE

- 1 Turn the button to the "pump" position, then depress the pump to inflate the chamber.
2. turn the button to the "let air out" position to deflate.

SIZE CHART (SHOE SIZE):

S: 36 - 39
M: 40 - 44
L: 44 - 46,5
XL: 46.5+

Length: AT53005-45 cm AT53006- 45 cm, AT53008 - 26 cm AT53030- 33 cm

COMPONENTS

Orthosis, instruction manual

NOTES:

In the event of a device-related "serious incident" that directly or indirectly led to, could have led to, or is likely to lead to any of the following:

- (a) death of a patient, user or other person, or
- (b) temporary or permanent deterioration of the health of a patient, user or other person, or
- (c) a serious risk to public health

the above 'serious incident' must be notified to the manufacturer and to the competent authority of the Member State in which the user or patient is resident. In the case of Poland, the competent authority is the Office for Registration of Medicinal Products, Medical Devices and Biocidal Products.

NOTES:

In case of pain, allergic reactions or other disturbing symptoms unclear to the user related to the use of the medical device, consult a healthcare professional.

NOTE: It is forbidden to use the product other than in accordance with its intended use

ATTENTION: The manufacturer is not responsible for damage caused by negligent maintenance, inadequate servicing or resulting from failure to follow the recommendations contained in this instruction manual.

CLEANING AND MAINTENANCE

Remove metal parts and casing before washing (clean the casing with a damp cloth and mild detergent). Wash by hand in cold water using mild soap. Dry flat only in the shade. Do not tumble dry, dry clean or use bleach. Do not expose to sources of heat.

STORAGE

Store the product in a dry cool room, away from moisture and direct sunlight.

HOW TO DISPOSE OF THE PRODUCT

The medical device can be disposed of as normal household waste.

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Model:
Serial number:

Point of sale stamp

and a legible signature of the Seller:

1. ANTAR Sp. J., 03-068 Warsaw, ul. Zawiślańska 43, hereby grants a 12-month warranty on the product, starting from the date of purchase by the Buyer.
2. During the warranty period ANTAR undertakes to carry out the repairs necessary to restore proper functioning of the product, within 14 days from the date of receipt of the product by the service centre, together with the proof of purchase and the warranty card.
3. ANTAR reserves the right to replace the product if the repair is deemed uneconomic. Repair or replacement of faulty products does not extend the warranty period.
4. If the product is found to be malfunctioning, please contact your retailer/distributor or (if it is impossible) send it to ANTAR (as specified in point 1).
5. Before use, read the instructions carefully and follow them. If the product is used improperly and contrary to recommendations, the warranty will be irrevocably lost.
6. The warranty does not cover damage caused by external mechanical force, load exceeding the maximum allowed weight of the user, as described in the manual, as well as contact of the product with alcohol, greasy substances or petrol. The warranty does not cover parts which are subject to normal wear or tear during use, such as the bearing, rubber parts and others.
7. The product delivered to the service in connection with the complaint cannot cause an epidemic hazard.
8. THIS WARRANTY IS ONLY VALID TOGETHER WITH THE PROOF OF PURCHASE (BILL, RECEIPT, INVOICE). IN CASE OF A COMPLAINT, A PROOF OF PURCHASE MUST BE ATTACHED.
9. Complaints submitted without the proof of purchase and a warranty card with the product serial number will not be considered.

ANTAR Sp. J.
03-068 Warsaw, ul. Zawiślańska 43
Tel. 22 518 36 00, Fax 22 518 36 30
www.antar.net

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