AT53069, AT53070, AT53071
KOŁNIERZ ORTOPEDYCZNY INSTRUKCJA UŻYTKOWANIA
ORTHOPEDIC COLLAR USER MANUAL
ORTHOPÄDISCHE HALSKRAUSE GEBRAUCHSANWEISUNG
ORTOPEDICKÝ LÍMEC NÁVOD K POUŽITÍ
ORTOPEDICKÝ GOLIER NÁVOD NA POUŽITIE
COLLIER ORTHOPÉDIQUE MANUEL D'UTILISATION
ORTHOPEDISCHE HALSKRAAG GEBRUIKSAANWIJZING
COLLAR ORTOPÉDICO MANUAL DE USO
COLLARE ORTOPEDICO MANUALE D'USO
ORTOPEDISK KRAGE BRUKSANVISNING

MD

PL

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EN

Thank you for purchasing our product. Before use, please carefully read this user manual. Improper usage may result in health risks, bodily injuries, or accidents. Therefore, to ensure the safe and effective use of the product, please thoroughly familiarize yourself with this manual.

WARNING

Check all parts for damage that may have occurred during transport. If any damage is found, DO NOT use the product. For further information, contact the manufacturer.

APPLICATION/INDICATIONS AT53069:

Neck stabilization
Neck sprains and strains
Mild cervical osteoarthritis
Acute or chronic painful cervical syndrome

FEATURES OF AT53069

The soft collar is made of low-density foam covered with an additional material for maximum patient comfort. A Velcro closure allows easy adjustment. Contoured foam provides comfortable, anatomical fitting. The rounded foam edge toward the head ensures additional wearing comfort.

APPLICATION/INDICATIONS AT53070:

Neck stabilization
Mild cervical osteoarthritis
Acute or chronic cervical spine pain syndrome
Postoperative treatment
Rheumatic diseases of the cervical spine

FEATURES OF AT53070

The semi-rigid collar is made of medium-density foam covered with additional material for maximum patient comfort. A Velcro closure allows easy adjustment. Contoured foam ensures a comfortable, anatomical fit with added support provided by plastic inserts. The rounded foam edge toward the head ensures additional user comfort.

USER MANUAL AT53069/AT53070

Applying the cervical collar

Positioning:

Place the collar around the neck with the rounded foam edge facing upwards to support the chin.
Ensure the collar fits evenly against the skin without leaving any gaps.

Fastening:

Fasten the collar using the Velcro closure at the back.
Adjust the closure to ensure the collar is stable but not too tight to compress the neck or restrict free breathing.

During use

Regular checks:

Ensure the collar does not cause discomfort, chafing, or pain.
Check that the collar is neither too loose nor too tight.

Medical recommendations:

Wear the collar according to the instructions of your doctor or physiotherapist, particularly regarding the duration of use.

Daily activities:

Avoid sudden movements of the neck or head to maintain stabilization.
If you experience discomfort or pain, consult a doctor.

Removing the cervical collar

Loosening the closure:

Undo the Velcro closure at the back of the collar.

Careful removal:

Gently remove the collar, avoiding sudden movements.

APPLICATION/INDICATIONS AT53071:

Cervical spondylosis
Postoperative cervical spine treatment
Cervical spine fractures
Rheumatic diseases of the cervical spine

FEATURES OF AT53071

The rigid collar is made of high-density foam. The lightweight, two-piece collar is equipped with rigid occipital and mandibular supports. Its front and rear attachments further restrict cervical movements. Additionally, the collar features a tracheal opening.

USER MANUAL AT53071

Preparation:

Ensure the collar is clean and dry.

Positioning:

Place the back part of the collar behind your neck and the front part against the front of your neck, ensuring the tracheal opening is positioned centrally under your chin.

Adjust the lower edge of the collar to align with your collarbones and the upper edge to align with the bottom of your jawline.

Fastening:

Connect the front and back parts using the Velcro fasteners on both sides of the neck.
Ensure the collar is stable but not too tight to compress the neck or restrict breathing.

During use:

Regularly check that the collar is well-fitted and does not cause chafing, pain, or discomfort.

Avoid sudden movements of the neck and head to maintain stabilization.
 Use the collar as recommended by your doctor, particularly concerning the duration of use.
 Removing the collar:

Unfasten the Velcro closures on both sides of the neck.
 Gently remove the collar, avoiding sudden movements.

CONTRAINDICATIONS

Physical or mental limitations (e.g., visual impairment) that prevent safe handling of the product.
 In cases where cervical spine fractures or injuries are severe and unstable, a simple orthopedic collar may not provide sufficient stabilization. Such cases require more advanced immobilization solutions or surgical intervention.
 The use of an orthopedic collar may worsen existing skin conditions on the neck, such as eczema, psoriasis, or infections. Lack of proper hygiene and limited ventilation under the collar can exacerbate skin problems.
 Caution is advised for individuals with circulatory disorders in the neck. Excessive collar tightness may impair blood flow and lead to further complications.
 The cervical collar restricts neck mobility, which may lead to muscle atrophy and sensory impairment. Regular physical evaluations are recommended, especially during prolonged use.
 Prolonged use without regular medical consultations may weaken neck muscles and worsen the condition. Always use the collar under medical supervision.

TARGET PATIENT GROUP

Individuals suffering from diseases, dysfunctions, or injuries requiring treatment/rehabilitation or compensation for which this product is intended (see the product's intended purpose in this manual). The product can be purchased independently or based on the recommendation of a physician, therapist, or other specialist. In both cases, factors such as available sizes, necessary features, dimensions, variations, indications, and contraindications must be considered, as well as the information provided by the manufacturer.

WARNING: Do not use the product after applying warming ointments or gels.

SIZES

To select the appropriate size, measure the neck circumference.

AT53069/AT53070:

S: 28–33 cm
 M: 34–39 cm
 L: 40–45 cm
 Height: 8 cm

AT53071:

XS: 28–32 cm
 S: 33–37 cm
 M: 38–42 cm
 L: 43–47 cm
 Height: 13 cm

CAUTION: Using the product in ways other than intended is prohibited.

NOTE:

In the event of a "serious incident" related to the product that directly or indirectly resulted, could have resulted, or may result in any of the following:

- a) death of a patient, user, or other individual;
- b) temporary or permanent deterioration of health of the patient, user, or another individual;
- c) serious public health threat;

such a "serious incident" must be reported to the manufacturer and the competent authority in the user's or patient's country of residence. In Poland, the competent authority is the Office for Registration of Medicinal Products, Medical Devices, and Biocidal Products.

WARNING:

If pain, allergic reactions, or other unclear or concerning symptoms occur during use, discontinue use immediately and consult a healthcare professional.

CLEANING INSTRUCTIONS:

Hand wash in lukewarm water with soap. Do not use bleach or harsh detergents. Do not iron, tumble dry, or dry clean. Store in a cool, dry place away from direct sunlight.

	Hand wash at up to 40 degrees C
	Do not bleach
	Do not iron
	Do not tumble dry
	Do not dry clean

DISPOSAL OF THE PRODUCT

Dispose of the product in accordance with applicable waste disposal regulations. Once withdrawn from use, the medical device can be disposed of as regular municipal waste.

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Model:
Serial number:

Point of sale stamp

and a legible signature of the Seller:

1. ANTAR Sp. J., 03-068 Warsaw, ul. Zawiślańska 43, hereby grants a 12-month warranty on the product, starting from the date of purchase by the Buyer.
2. During the warranty period ANTAR undertakes to carry out the repairs necessary to restore proper functioning of the product, within 14 days from the date of receipt of the product by the service centre, together with the proof of purchase and the warranty card.
3. ANTAR reserves the right to replace the product if the repair is deemed uneconomic. Repair or replacement of faulty products does not extend the warranty period.
4. If the product is found to be malfunctioning, please contact your retailer/distributor or (if it is impossible) send it to ANTAR (as specified in point 1).
5. Before use, read the instructions carefully and follow them. If the product is used improperly and contrary to recommendations, the warranty will be irrevocably lost.
6. The warranty does not cover damage caused by external mechanical force, load exceeding the maximum allowed weight of the user, as described in the manual, as well as contact of the product with alcohol, greasy substances or petrol. The warranty does not cover parts which are subject to normal wear or tear during use, such as the bearing, rubber parts and others.
7. The product delivered to the service in connection with the complaint cannot cause an epidemic hazard.
8. THIS WARRANTY IS ONLY VALID TOGETHER WITH THE PROOF OF PURCHASE (BILL, RECEIPT, INVOICE). IN CASE OF A COMPLAINT, A PROOF OF PURCHASE MUST BE ATTACHED.
9. Complaints submitted without the proof of purchase and a warranty card with the product serial number will not be considered.

ANTAR Sp. J.
03-068 Warsaw, ul. Zawiślańska 43
Tel. 22 518 36 00, Fax 22 518 36 30
www.antar.net

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