
AT92022

**INSTRUKCJA OBSŁUGI ORTEZA KOLANA
KNEE ORTHOSIS USER MANUAL
KNIEORTHESE GEBRAUCHSANLEITUNG
NÁVOD K POUŽITÍ KOLENNÍ ORTÉZY
NÁVOD NA POUŽITIE ORTÉZA KOLENA
MANUEL D'UTILISATION DE L'ORTHÈSE DE GENOU
KNIE-ORTHOSE GEBRUIKERSHANDLEIDING
MANUAL DEL USUARIO DE LA ORTESIS DE RODILLA
MANUALE D'USO DELL'ORTOSI DEL GINOCCHIO
KNÄORTOS ANVÄNDARHANDBOK**

PL

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ENG

PURPOSE / INDICATIONS

The brace is made of a 4-way stretchable material. It strengthens the knee joint, acts prophylactically, protecting it against injuries, especially when the joint is weakened, in the case of rheumatic and degenerative diseases.

CONTRAINDICATIONS

physical or mental limitations (e.g. visual impairment) that prevent the safe handling of the product.

WARNING

This product is for personal use only. Do not use the product if there are any damage to the included accessories, components or the product itself. Do not wear the orthosis if ointments or warming lotions have been applied to the area. If a rash or irritation occurs, remove the orthosis and discontinue use, consult a physician immediately. The product should be used according to the doctor's instructions.

INSTRUCTIONS FOR USE

The product may be available in various sizes. Please refer to the size chart to select the correct size.

WASHING INSTRUCTIONS

Hand wash with mild soap at 40 ° C, Air dry away from heat and direct sunlight, Do not bleach, iron or dry clean.

STORAGE

Store in a dry place away from the sun.

SIZES (measure the circumference of the knee on a straight leg):

S: 28,6 - 33,7 cm

M: 33,7 - 40,6 cm

L: 40,6 - 47,6 cm

XL: 47,6 - 55,2 cm

XXL: 55,2 - 62,9 cm

CAUTION:

When a product-related "serious incident" occurs, which directly or indirectly led, could or may have led to any of the following events:

- death of a patient, user or other person or
- temporary or permanent deterioration in the health of a patient, user or other person, or
- a serious threat to public health

this "serious incident" must be reported to the manufacturer and to the competent authority of the Member State in which the user or patient is resident. In the case of Poland, the competent authority is the Office for Registration of Medicinal Products, Medical Devices and Biocidal Products.

CAUTION:

In the event of pain, allergic reactions or other worrying symptoms, unclear for the user, related to the use of a medical device, please consult a healthcare professional.

ATTENTION: it is forbidden to use the product other than in accordance with its intended purpose

ATTENTION: The manufacturer is not responsible for damage caused by neglect of maintenance, inadequate servicing or resulting from failure to comply with the recommendations contained in this manual.

METHOD OF DISPOSAL OF THE PRODUCT AFTER REMOVING FROM SERVICE

After the device is withdrawn from use, the medical device can be disposed of as normal household waste.

SCOPE OF DELIVERY:

orthosis, manual.

TARGET GROUP OF PATIENTS

The healthcare professional shall indicate the use of the device to adults and children on his/her own responsibility, taking into account the available variants/sizes/necessary functions/sizes and indications, bearing in mind the information provided by the manufacturer.

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WARRANTY CARD

Model:
Serial number:
Point of sale stamp
and a legible signature of the Seller:

1. ANTAR Sp. J., 03-068 Warsaw, ul. Zawisłańska 43, hereby grants a 12-month warranty on the roller, starting from the date of purchase by the Buyer. Products financed by the National Health Fund are covered with the extended warranty period, which covers half of the useful life specified in the Regulation on Guaranteed Services, provided that the useful life is specified in the Regulation in months or years.
2. During the warranty period ANTAR undertakes to carry out the repairs necessary to restore proper functioning of the product, within 14 days from the date of receipt of the product by the service centre, together with the proof of purchase and the warranty card.
3. ANTAR reserves the right to replace the product if the repair is deemed uneconomic. Repair or replacement of faulty products does not extend the warranty period.
4. If the product is found to be malfunctioning, send it immediately to ANTAR (as specified in point 1) or contact the point of sale.
5. Before use, read the instructions carefully and follow them. If the product is used improperly and contrary to recommendations, the warranty will be irrevocably lost.
6. The warranty does not cover damage caused by external mechanical force, load exceeding the maximum permissible weight of the user as described in the manual, as well as contact of the product with alcohol, greasy substances or petrol. The warranty does not cover parts subject to normal wear and tear during use, such as the bearing, rubber parts and others.
7. The product delivered to the service in connection with the complaint cannot cause an epidemic.
8. THIS WARRANTY IS ONLY VALID TOGETHER WITH THE PROOF OF PURCHASE (BILL, RECEIPT, INVOICE). IN CASE OF A COMPLAINT, A PROOF OF PURCHASE MUST BE ATTACHED. IN ORDER TO EXTEND WARRANTY PERIOD FOR PRODUCTS FINANCED BY THE NATIONAL HEALTH FUND, ATTACH A COPY OF THE ORDER FOR THE SUPPLY OF MEDICAL DEVICES BEING ORTHOPAEDIC ITEMS AND AIDS.
9. Complaints submitted without the proof of purchase and a warranty card with the product serial number will not be considered.

ANTAR Sp. J.
03-068 Warsaw, ul. Zawisłańska 43
Tel. 22 518 36 00, Fax 22 518 36 30
www.antar.net

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