
AT92068, AT92164

**INSTRUKCJA OBSŁUGI ORTEZA KRĘGOSŁUPA
USER MANUAL SPINAL ORTHOSIS
GEBRAUCHSANLEITUNG WIRBELSÄULENORTHESSE
UŽIVATELSKÁ PŘÍRUČKA BEDERNÍ ORTÉZA
NÁVOD NA POUŽITIE ORTÉZA CHRBTICE
MANUEL D'UTILISATION ORTHÈSE RACHIDIENNE
GEBRUIKERSHANDLEIDING SPINALE ORTHOSE
MANUAL DEL USUARIO ORTESIS ESPINAL
MANUALE D'USO ORTOSI SPINALE
ANVÄNDARINSTRUKTIONER RYGGGORTOS**

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PURPOSE / INDICATIONS

Product is used to strengthen the lumbar spine. Sewn-in reinforcements, as well as side supports provide additional stiffening. Flexible, "breathable" material allows you to maintain natural body temperature and improves blood circulation. The belt is recommended for various diseases of the spine, e.g. lumbar discopathy; osteoarthritis of the spine; root syndromes (sciatica).

CONTRAINDICATIONS

physical or mental limitations (e.g. visual impairment) that prevent the safe handling of the product.

WARNING

This product is for personal use only. Do not use the product if there are any damages on the attached accessories, components, or at the product itself. Do not wear the orthosis if ointments or warming lotions have been applied to the area. If a rash or irritation occurs, remove the orthosis and discontinue use. Consult a physician immediately. The product should be used according to the doctor's instructions.

INSTRUCTIONS FOR USE

The product may be available in various sizes. Please refer to the size chart to select the correct size. Take the product out of the packaging and then put it on the appropriate part of the body.

WASHING INSTRUCTIONS

Hand wash with mild soap at 40 ° C. Air dry away from heat and direct sunlight. Do not bleach, iron or dry clean. Before washing, remove the stiffeners.

STORAGE

Store in a dry place away from the sun.

SIZE AT92068 (waist circumference):

S: 58.4 - 71.1 cm
 M: 71.1 - 83.8 cm
 L: 83.8 - 96.5 cm
 XL: 96.5 - 109.2 cm
 XXL: 109.2 - 121.9 cm
 XXXL: 121.9 - 134.6 cm

SIZE AT92164 (waist circumference):

S: 58.4 - 71.1 cm
 M: 71.1 - 83.8 cm
 L: 83.8 - 96.5 cm
 XL: 96.5 - 109.2 cm
 XXL: 109.2 - 121.9 cm
 XXXL: 121.9 - 139.7 cm

CAUTION:

When a product-related "serious incident" occurs, which directly or indirectly led, could or may have led to any of the following events:

- death of a patient, user or other person or
- temporary or permanent deterioration in the health of a patient, user or other person, or
- a serious threat to public health

this "serious incident" must be reported to the manufacturer and to the competent authority of the Member State in which the user or patient is resident. In case of Poland, the competent authority is the Office for Registration of Medicinal Products, Medical Devices and Biocidal Products.

CAUTION:

In the event of pain, allergic reactions or other disturbing, unclear for the user symptoms related to the use of a medical device, please consult a healthcare professional.

ATTENTION: it is forbidden to use the product other than in accordance with its intended purpose

HOW TO DISPOSE OF THE PRODUCT

The medical device can be disposed of as normal municipal waste.

SCOPE OF DELIVERY: orthosis, manual

TARGET GROUP OF PATIENTS

Persons suffering from diseases, dysfunctions or injuries for the treatment/rehabilitation or compensation of which this device is intended (see section on intended use of the device in these instructions). The device can be purchased by the user on their own or on the recommendation of a doctor, therapist or other specialist. Whether you purchase the device yourself or on the recommendation of a doctor/therapist/other specialist, you must take into account the available sizes/necessary functions/size/variants of the device, the indications and contraindications for use, as well as the information provided by the manufacturer.

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Model:

Serial number:

Point of sale stamp
and a legible signature of the Seller:

1. AN TAR Sp. J., 03-068 Warsaw, ul. Zawiślńska 43, hereby grants a 12-month warranty on the product, starting from the date of purchase by the Buyer.
2. During the warranty period AN TAR undertakes to carry out the repairs necessary to restore proper functioning of the product, within 14 days from the date of receipt of the product by the service centre, together with the proof of purchase and the warranty card.
3. AN TAR reserves the right to replace the product if the repair is deemed uneconomic. Repair or replacement of faulty products does not extend the warranty period.
4. If the product is found to be malfunctioning, please contact your retailer/distributor or (if it is impossible) send it to AN TAR (as specified in point 1).
5. Before use, read the instructions carefully and follow them. If the product is used improperly and contrary to recommendations, the warranty will be irrevocably lost.
6. The warranty does not cover damage caused by external mechanical force, load exceeding the maximum allowed weight of the user, as described in the manual, as well as contact of the product with alcohol, greasy substances or petrol. The warranty does not cover parts which are subject to normal wear or tear during use, such as the bearing, rubber parts and others.
7. The product delivered to the service in connection with the complaint cannot cause an epidemic hazard.
8. THIS WARRANTY IS ONLY VALID TOGETHER WITH THE PROOF OF PURCHASE (BILL, RECEIPT, INVOICE). IN CASE OF A COMPLAINT, A PROOF OF PURCHASE MUST BE ATTACHED.
9. Complaints submitted without the proof of purchase and a warranty card with the product serial number will not be considered.

ANTAR Sp. J.
03-068 Warsaw, ul. Zawiślńska 43
Tel. 22 518 36 00, Fax 22 518 36 30
www.antar.net

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