AT92204
INSTRUKCJA UŻYTKOWANIA, ORTEZA STAWU SKOKOWEGO
INSTRUCTION MANUAL, ANKLE ORTHOSIS
GEBRAUCHSANWEISUNG, KNÖCHELSCHIENE
NÁVOD K POUŽITÍ, KOTNÍKOVÁ ORTÉZA
NÁVOD NA POUŽITIE, ORTÉZA NA ČLENOK
MANUEL D'INSTRUCTION, ORTHÈSE DE CHEVILLE
GEBRUIKSAANWIJZING, ENKELBRACE
MANUAL DE INSTRUCCIONES, TOBILLERA
MANUALE DI ISTRUZIONI, CAVIGLIERA
BRUKSANVISNING, FOTLEDSSTÖD

MD

CE

PL

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EN

Thank you for purchasing our product, please read the instructions for use carefully. Do not use the product before reading and understanding this manual. If you do not understand the warnings, cautions and recommendations, please contact your health care professional or dealer to avoid product damage or personal injury.

CAUTION: Inspect all parts for damage caused in transit. If you notice such damage, DO NOT use the product. Contact the manufacturer for more information.

INDICATIONS

The orthosis is made of flexible, breathable material. It strengthens the ankle joint, acts prophylactically to protect it from injury, especially when the joint is weakened, in case of rheumatic and degenerative diseases. It is also helpful during post-injury rehabilitation.

CONTRAINDICATIONS

Physical or mental limitations (such as visual impairment) that prevent safe handling.

- open wounds
- severe swelling and deformities

WARNING

This product is intended for personal use only. Do not use the product if there is any damage on the supplied accessories, components or on the product itself. Do not wear the orthosis if ointments or warming fluids have been applied to the area. If a rash or irritation develops, remove the orthosis and discontinue use, consult a doctor immediately. The product should be used as directed by your doctor.

INSTRUCTIONS FOR USE

The product may be available in different sizes. Refer to the size chart to select the appropriate size. The product should be slipped over the foot, securing the ankle joint properly.

WASHING INSTRUCTIONS

Hand wash with mild soap at up to 40 °C. Air dry laid flat away from heat and direct sunlight. Do not use strong detergents. Do not bleach. Do not iron. Do not dry clean.

STORAGE

Store in a dry and cool place, away from direct sunlight.

SIZE (circumference around the ankle joint):

- S 15.6cm - 21.0cm
- M 21,6cm - 26,0cm
- L 26.7cm - 31.1cm
- XL 31,8cm - 36,3cm
- XXL 36.8cm - 41.3cm

NOTE

In the event of a product-related "serious incident" that directly or indirectly led to, may have led to, or may lead to any of the following:

- (a) death of a patient, user or other person, or
- (b) temporary or permanent deterioration of the health of a patient, user or other person, or

(c) a serious threat to public health
the above "serious incident" must be reported to the manufacturer and to the competent authority of the Member State in which the user or patient resides. In the case of Poland, the competent authority is the Office for Registration of Medicinal Products Medicinal Products, Medical Devices and Biocidal Products.

NOTE

In case of pain, allergic reactions or other disturbing symptoms unclear to the user related to the use of the medical device, consult a health care professional.

CAUTION

It is forbidden to use the product other than in accordance with its intended use

NOTE

The manufacturer is not responsible for damages caused by negligent maintenance, improper service or resulting from failure to follow the recommendations contained in this instruction manual.

METHOD OF DISPOSAL OF THE PRODUCT AFTER ITS WITHDRAWAL FROM SERVICE

After the device is decommissioned, the medical device may be disposed of as ordinary municipal waste, with the exception of electrical products - proceed in a manner appropriate for the disposal of electrical and electronic equipment.

COMPONENTS

Orthosis, instructions for use.

TARGET PATIENT GROUP

Persons struggling with diseases, dysfunctions or injuries, for the treatment/rehabilitation or compensation of which this product is intended (see the section on the purpose of the product of this manual). The device may be purchased by the user independently or on the recommendation of a physician, therapist or other specialist. Whether purchasing the device on its own or as indicated by a doctor/therapist/other specialist, the available sizes/necessary functions and variants of the device, indications and contraindications for use, as well as information provided by the manufacturer, should be taken into account.

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Model:
Serial number:

Point of sale stamp

and a legible signature of the Seller:

1. ANTAR Sp. J., 03-068 Warsaw, ul. Zawisłańska 43, hereby grants a 12-month warranty on the product, starting from the date of purchase by the Buyer.
2. During the warranty period ANTAR undertakes to carry out the repairs necessary to restore proper functioning of the product, within 14 days from the date of receipt of the product by the service centre, together with the proof of purchase and the warranty card.
3. ANTAR reserves the right to replace the product if the repair is deemed uneconomic. Repair or replacement of faulty products does not extend the warranty period.
4. If the product is found to be malfunctioning, please contact your retailer/distributor or (if it is impossible) send it to ANTAR (as specified in point 1).
5. Before use, read the instructions carefully and follow them. If the product is used improperly and contrary to recommendations, the warranty will be irrevocably lost.
6. The warranty does not cover damage caused by external mechanical force, load exceeding the maximum allowed weight of the user, as described in the manual, as well as contact of the product with alcohol, greasy substances or petrol. The warranty does not cover parts which are subject to normal wear or tear during use, such as the bearing, rubber parts and others.
7. The product delivered to the service in connection with the complaint cannot cause an epidemic hazard.
8. THIS WARRANTY IS ONLY VALID TOGETHER WITH THE PROOF OF PURCHASE (BILL, RECEIPT, INVOICE). IN CASE OF A COMPLAINT, A PROOF OF PURCHASE MUST BE ATTACHED.
9. Complaints submitted without the proof of purchase and a warranty card with the product serial number will not be considered.

ANTAR Sp. J.
03-068 Warsaw, ul. Zawisłańska 43
Tel. 22 518 36 00, Fax 22 518 36 30
www.antar.net

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