HF6001, HF6002, HF6008

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Thank you for purchasing the mattresses we offer. Please read these instructions before use.

To avoid an accident please:

1. connect the pump to the correct type of electrical outlet (220-230V / 50 Hz).
2. remove the plug from the socket whenever the pump is not running.
3. do not position the pump where it could come into contact with water.
4. not touching the pump and disconnecting it from the mains if it gets wet.
- (5) Please keep your proof of purchase of the mattress (receipt, bill or invoice). You will need it in case you need to contact our service.

To avoid injury and fire:

1. the pump should never be left running unattended.
2. constant supervision is necessary if the pump is used near children or invalids.
- 3 The pump may only be used for the purpose for which it is intended. Do not attach unauthorised accessories.
- (4) Do not use a pump with a damaged cord or plug if it is not working properly, has been dropped from a great height, damaged or wet. Such a pump must be returned to the point of purchase, who will forward it for service.
5. keep the cable away from hot surfaces.
6. do not position the pump where the cable or hoses/tubes may present a hazard to persons unfamiliar with their use.

When disconnecting from a socket, please pull the plug, not the cable.

8 Do not clog the pump filter. Do not place the pump on a soft surface where the filter could become blocked.

The filter must be free of debris such as hair, dander etc.

9. consult your doctor before use.

Preparation for use and use

Note: The mattress is fitted with a pad specially adapted to the pump. Use of any other pad may result in malfunction of the set.

- 1 Place the mattress on the bed. Pay attention to the direction mattress. The air tube inlet points should be located is at your feet.
2. check that all chambers are connected (does not apply to bubble model).
- 3 Set the switch to the off - "0" position.
4. connect the pump to the mains.
5. connect the tubing to the pump by pushing it firmly into place into the air outlets located on the side of the pump.
- 6 Place the pump on a firm surface or hang it from the hinges on the bottom of the pump.
- 7 The location where the pump is set up must be clean. Dirt, such as dust or hair, can block the filter.
- 8 Arrange the tubes so that they are not pinched or kinked.
9. turn the switch to the on position - "1".
10. Set the pressure regulator to maximum. This will reduce the pumping time of the mattress, which takes about 15- 20 min.
11. Inflate the mattress to the doctor-recommended pressure or the most comfortable level for the patient. Ensure the mattress is not too soft when pressure is lowered, and prevent the patient from pressing body parts against the mattress base. If unsure about the appropriate pressure, leaving the mattress at the maximum pressure is an option.
12. Some tube mattress models feature a "STATIC" function. In this mode, all chambers receive constant air pressure, deactivating alternating pressure mode. Use the "STATIC" function for specific purposes (e.g., washing or transporting patients), but note that it reduces the mattress's anti-decubitus properties.

NOTE: When connecting the pump to a mattress with STATIC function, leave the centre hole in the pump free.

NOTE: If the patient is on the mattress and no activities are carried out around the patient, the „STATIC" function must be OFF. Failure to do so may result in pressure sores.

Consult your doctor for precise usage instructions and optimal pressure. The set induces a 'continuous motion' effect, promoting pressure changes on specific body areas, continuous muscle massage, and improved blood flow. This mattress serves as a valuable tool for pressure sore prevention. Typically, alternate chambers inflate and empty, with a reversal occurring after about 5 - 10 minutes.

**NOTE: DO NOT MOVE THE MATTRESS (INFLATED/UNINFLATED) WHILE THE PATIENT IS ON IT.
ATTENTION: THIRD-PERSONS SHOULD NOT BE SITTING (sitting, lying) ON THE INFLATED MATERIAL WHEN A SICK PERSON IS ON IT.**

NOTE: DURING THE ENTIRE PERIOD WHEN THE PATIENT IS ON THE MATTRESS, THE PUMP MUST BE IN OPERATION!

This is a condition of alternating filling and emptying of the chambers every 5-10 minutes or so.

CAUTION: SWITCHING OFF THE PUMP WHEN THE PATIENT IS ON THE MATTRESS WILL QUICKLY LEAD TO PRESSURE SORES AND ENDANGER THE PATIENT'S HEALTH!

NOTE: IF YOU FIND THAT THE CHAMBERS ARE NOT ALTERNATELY FILLED WITH AIR AND EMPTIED EVERY APPROX. 5-10 MINUTES, CONTACT THE SERVICE DEPARTMENT!

End of use (when the sick person is not on the mattress)

1. flip the switch to the off position.
2. disconnect the tubes from the air outlet to allow the remaining air in the mattress to escape.

Cleaning and storage

Disconnect the pump from the mains. Wipe all accessible surfaces with a wet cloth and liquid for cleaning and then with a cloth soaked in clean water. The pump must be completely dry before switching on.

Notes:

1. use of isopropyl alcohol or solvents may damage the pump.
2. Ensure the pump is dry before restarting. Wash the waterproof mattress cover (selected models) by hand or in a machine at max. 40°C, letting it dry completely. Clean the mattress (mattress pad) with soapy water or detergent, wiping it thoroughly and allowing complete drying. Do not boil, tumble-dry, bleach, dry-clean, iron, or compress the mattress. Avoid contact with alcohol or oily substances to prevent damage. UV light can be used for decontamination. Before storing, fold the mattress, ensuring it's entirely dry to prevent mold. Store the mattress in a dry, cool, shaded place.

Chamber replacement (for tube mattresses)

If one of the air chambers is damaged, it should be replaced with the spare included in the package.

1. To do this, switch off the pump, disconnect the cable from the mains socket.
2. let the air out of the mattress.
3. Detach the damaged air chamber.
4. Pull the air line out of the air line connector on either side of the damaged air chamber.
5. Remove the damaged air chamber from the mattress straps and replace it with a spare.
6. the air duct connector should be placed in the air ducts so that they fit properly, as before the ducts were disconnected.
7. Fasten the fitted air chamber.

Problem solving

If problems occur, please consult this section. If this does not help please contact your point of sale. Please do not attempt to repair the pump yourself or open it. There is a risk of electric shock.

I. The pump does not operate even though the switch is in the "on" position.

- a) Is the cable connected to the mains and also to the pump? If not, it should be connected.
- b) Is there voltage on the mains? If not, determine the reason for the lack of voltage.
- c) Is the cable damaged? If so, replace the cable.

II. The pump is working, but the mattress pumps slowly or does not reach the desired pressure.

- a) Is the pressure regulator set at the correct pressure? If not, it should be adjusted accordingly.
- b) Are the tubes inserted properly into the air outlets on the pump? If not, tighten them so that there is no air loss.
- c) Are the tubes bent or crushed? If so, straighten the tubes.
- d) Is one of the chambers leaking air? If so, depending on the model, seal the hole with the supplied patch, replace the chamber with a spare, and contact the dealer if the spare parts are not available (out of stock).
- e) Is there a loss of air inside the pump (loud hissing sound)? If so, return the pump for repair.

CONTRAINDICATIONS

Physical or mental limitations (e.g. visual impairment) that prevent safe handling . Contraindications to remaining in a supine position.

TARGET PATIENT GROUP

Individuals with conditions for which this device is intended (see the instructions' intended use section) can buy it independently or on professional recommendation. Whether self-purchased or recommended by a doctor/ therapist/specialist, consider device sizes, functions, variants, indications, contraindications, and manufacturer-provided information.

NOTES:

In the event of a product-related "serious incident" that directly or indirectly led to, could have led to or may lead to any of the following:

- (a) the death of a patient, user or other person; or
- (b) temporary or permanent deterioration of the health of a patient, user or other person, or
- (c) a serious threat to public health

Report the mentioned 'serious incident' to the manufacturer and the competent authority in the user or patient's resident Member State. In Poland, contact the Office for Registration of Medicinal Products, Medical Devices, and Biocidal Products.

NOTES:

If you experience pain, allergic reactions or other distressing, unclear symptoms related to the use of the medical device, stop using the device and consult a healthcare professional.

HOW TO DISPOSE OF THE PRODUCT

This medical device is an electrical product - handle in a manner appropriate for the disposal of electrical and electronic equipment.

NOTE: Before use, check that all components have been assembled correctly. It is forbidden to use the product other than for its intended purpose.

(HF6002/HF6P01) E TUBULAR MATTRESS XI (ALU 18 R)

Pneumatic tubular mattress with aluminum pump and pressure control, offering effective pressure sore prevention for bedridden, long-term, or terminally immobile patients, including those recovering from strokes, paralyzed individuals, ICU patients, and neurological cases.

Features of the tube mattress:

- detachable tubes, ventilation function, cycle - 12 minutes.
- CPR valve for rapid deflation of the mattress in case of emergency
- maximum user weight: 150 kg

Primer:

- dimensions: 200 x 90 x 11.5 cm - 18 cells
- composition of the mattress cover: PU cover
- Air chamber composition: Nylon PVC
- repair kit and hose included, with additional chamber
- mattress cover
- underlay weight: 5.5 -6.5 kg

Pump (aluminium):

- power supply: 220V AC / 50Hz / 10W
- capacity: 6-8 litres/min.
- pressure range: 50 mmHg - 120 mmHg
- cycle time: 12 min
- size: 27.5 x 12 x 10 cm
- weight: 1.1 kg
- noise level 30 dB

(HF6002-17P/HF6P01) E TUBULAR MATTRESS XI (ALU)

Pneumatic tubular mattress with aluminum pump and pressure control, offering effective pressure sore prevention for bedridden, long-term, or terminally immobile patients, such as those recovering from strokes, paralyzed individuals, ICU patients, and neurological cases.

Features of the tube mattress:

- detachable tubes, ventilation function, cycle - 8 minutes
- CPR valve for quick deflation of the mattress in case of emergency
- maximum user weight: 150 kg

Primer:

- dimensions: 190 x 85 x 11.5 cm - 17 cells
- 4 ventilated tubes
- composition of mattress cover: nylon PU
- Air chamber composition: PVC
- repair kit and hose included, with additional chamber
- mattress cover
- backing weight: 5-5.5 kg

Pump (aluminium):

- power supply: 220V AC / 50Hz / 10W
- capacity: 6-8 litres/min.
- pressure range: 50 mmHg - 120 mmHg
- cycle time: 12 min
- size: 27.5 x 12 x 10 cm
- weight: 1.1 kg
- noise level 30 dB

(HF6002-22P/HFP01ABS) E TUBULAR MATTRESS XI (EKO)

Pneumatic mattress with alternating pressure, tubular design, and aluminum pump with pressure control. Offers effective protection against pressure sores for bedridden, chronically, or terminally immobile individuals, such as those recovering from strokes, paralyzed individuals, ICU patients, and neurological cases.

Features of the tube mattress:

- detachable tubes, cycle - 12 minutes.
- maximum user weight: 150 kg

Primer:

- dimensions: 186 x 80 x 8.5 cm - 22 cells
- composition of air chamber: PVC
- backing weight: 4.5-5 kg
- repair kit and hose included

Pump (plastic)

- power supply: 220V AC / 50Hz / 10W
- capacity: 5-7 litres/min.
- pressure range: 50 mmHg - 120 mmHg
- cycle time: 12 min
- size: 27.5 x 12 x 10 cm
- weight: 1 kg
- noise level 40 dB

(HF6008) Alternating pressure mattress with pump

Alternating pressure mattress with pump, tubular. Provides adequate protection for bedridden patients, long term or temporarily immobile, e.g. paralysed persons, patients after a stroke, neurological patients and people in the ICU.

Features of the tube mattress:

detachable tubes, ventilation function, cycle - 12 minutes
CPR valve for rapid deflation of the mattress in case of emergency
maximum user weight: 165 kg

Primer:

- dimensions: 200 x 90 x 20.5 - 20 cells
- 8 ventilated tubes
- composition of mattress cover: nylon PU
- Air chamber composition: Nylon TPU
- repair kit and hose included, with additional chamber
- mattress cover
- backing weight: 8 kg

Pump (aluminium):

- power supply: 220V/50Hz
- capacity: 6-8 litres / min
- pressure range: 50 mmHg - 120 mmHg
- cycle time: 12 min
- size: 27.5 x 13.5 x 9 cm
- weight: 1.6 kg
- noise level 30 dB

The pump has a CPC, it can retain air when the power is switched off. The pump has a low pressure alarm system.

BUBBLE MATTRESS XI(N)

Inflatable, alternating-pressure, bubble mattress, pressure-controlled pump. Comfortable and economical A mattress that provides good protection against the development of pressure sores. Minimises vibrations.

Primer:

- Dimensions (mattress inflated): 200 cm (L) x 90 cm (W) x 7 cm (H)
- material: medical grade PVC, highly elastic
- thickness of PVC: 0.30 mm
- repair kit and hose included
- 1.5 m hose
- 130 chambers
- weight: 1.8 kg
- maximum user weight: 130 kg

Pump (plastic)

- dimensions: 27.5 x 12 x 10 cm
- power supply: 220V/50Hz, 7W
- capacity: 5-7 l/min.
- 50 mmHg - 110 mmHg
- cycle time: 12 min.
- weight: 1 kg
- noise level 40 dB

BUBBLE MATTRESS XI(ALU)

Inflatable, alternating-pressure, bubble mattress, pressure-controlled pump. Comfortable and economical. Mattress provides good protection against pressure sores. Extremely quiet, aluminium pump motor cover. Minimising vibration.

Primer:

- Dimensions (mattress inflated): 200 cm (L) x 90 cm (W) x 7 cm (H)
- material: medical-grade PVC, highly flexible
- thickness of PVC: 0.30 mm
- repair kit and hose included
- 1.5 m hose
- 130 chambers
- weight: 1.8 kg
- maximum user weight: 130 kg

Pump (aluminium insert, very quiet):

- dimensions: 27.5 x 12 x 10 cm
- power supply: 220V/50Hz, 7W
- capacity: 6-8 l/min.
- 50mmHg - 110 mmHg
- cycle time: 12 min.
- weight: 1.1 kg
- noise level 30 dB

THIS SIGN INDICATES THE MAXIMUM WEIGHT OF THE USER.

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WARRANTY CARD

Model:
Serial number:

Point of sale stamp
and a legible signature of the Seller:

1. ANTAR Sp. J., 03-068 Warsaw, ul. Zawisłańska 43, hereby grants a 12-month warranty on the product, starting from the date of purchase by the Buyer.
2. During the warranty period ANTAR undertakes to carry out the repairs necessary to restore proper functioning of the product, within 14 days from the date of receipt of the product by the service centre, together with the proof of purchase and the warranty card.
3. ANTAR reserves the right to replace the product if the repair is deemed uneconomic. Repair or replacement of faulty products does not extend the warranty period.
4. If the product is found to be malfunctioning, please contact your retailer/distributor or (if it is impossible) send it to ANTAR (as specified in point 1).
5. Before use, read the instructions carefully and follow them. If the product is used improperly and contrary to recommendations, the warranty will be irrevocably lost.
6. The warranty does not cover damage caused by external mechanical force, load exceeding the maximum allowed weight of the user, as described in the manual, as well as contact of the product with alcohol, greasy substances or petrol. The warranty does not cover parts which are subject to normal wear or tear during use, such as the bearing, rubber parts and others.
7. The product delivered to the service in connection with the complaint cannot cause an epidemic hazard.
8. THIS WARRANTY IS ONLY VALID TOGETHER WITH THE PROOF OF PURCHASE (BILL, RECEIPT, INVOICE). IN CASE OF A COMPLAINT, A PROOF OF PURCHASE MUST BE ATTACHED.
9. Complaints submitted without the proof of purchase and a warranty card with the product serial number will not be considered.

FOSHAN HONGFENG CO.,LTD
Add. No.4-2 Leqiang Road,
Leping Sanshui Foshan,
528100, Guangdong, Chins

Authorized Representative:
MedNet EC-REP GmbH
Borkstasse 10 48163
Muenster Germany

IMPORTER:
ANTAR Sp. J. 03-068 Warsaw
ul. Zawisłanska 43 Tel. 22 518 36 00,
Fax 22 518 36 30 www.antar.net

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